PRODUCT CATALOG

DNA DIAGNOSTICS

Bacterial Infections
Viral Infections
Human Genetics

Only experts work with our **EliGene®** kits...

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Introduction

Dear customers,

ELISABETH PHARMACON is an innovative company developing advanced products for clinical and scientific laboratories. In our portfolio you will find plenty of unique polymerases, kits for DNA and RNA isolation, DNA diagnostics and human genetics that are presented in this catalog.

EliGene® kits satisfy only real experts

In the catalog we present a wide range of DNA/RNA diagnostic kits designed for bacterial and viral detection, as well as kits for human genetics and oncogenetics. The kits have been developed for laboratories that want to provide accurate results and will not be satisfied with lower accuracy of diagnostics.

EliGene® kits are universal

Design of our kits allows compatibility with a wide range of standard Real-Time PCR Cyclers from various manufacturers: ROCHE, Life Technologies, QIAGEN, etc.

EliGene® kits are exceptionally sensitive

High analytical sensitivity ensures laboratories to provide high quality diagnostics and doctors to have the right information to set up effective treatment. There is no false negativity with our kits, which may occur with some multiplexes. EliGene® kits detect one to ten copies of pathogenic DNA/RNA in the amplified sample.

EliGene® kits can provide quantitative results

Most of EliGene® DNA/RNA diagnostic kits are quantitative. They provide beside the information about presence of detected pathogens also information about their quantity.

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Content

Obsah

Kits for Bacterial Infections (Real-Time PCR kits)	3
Kits for Viral Infections (Real-Time PCR kits)	15
Kits for Human Genetics (Real-Time PCR kits)	27
Kits for Oncogenetics - (NGS)	31

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Molecular Beacons

Performance study and results

Within the frame of testing the functional characteristics of EliGene® Legionella pneumophila UNI, totally 50 clinical specimens were analyzed. The EliGene® Legionella pneumophila UNI Kit diagnosed correctly as *Legionella pneumophila* positive 5 specimens. Totally 45 specimens were rightly determined by EliGene® Legionella pneumophila UNI Kit as negative.

Sensitivity: 100%

Specificity: 100%

References / Zdroje

Heuner K, Swanson M. 2008. *Legionella: Molecular Microbiology*. Caister Academic Press. ISBN 978-1-904455-26-4.

Newton HJ, Ang DK, van Driel IR, Hartland EL. 2010. Molecular pathogenesis of infections caused by *Legionella pneumophila*. *Clin Microbiol Rev*. 23(2):274-98

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Molecular Beacons

Performance study and results

Within the frame of testing the functional characteristics of EliGene[®] Chlamydophila pneumoniae UNI Kit, totally 50 clinical specimens were analyzed. The EliGene[®] Chlamydophila pneumoniae UNI Kit diagnosed correctly as *Chlamydophila pneumoniae* positive 6 specimens. Totally 44 specimens were rightly determined by EliGene[®] Chlamydophila pneumoniae UNI Kit as negative.

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Specificita: 100 %

References / Zdroje

Pignanelli S, Shurdhi A, Delucca F, et al. 2009. Simultaneous use of direct and indirect diagnostic techniques in atypical respiratory infections from *Chlamydophila pneumoniae* and *Mycoplasma pneumoniae*. J Clin Lab Anal. 23(4):206–9.

Roulis E, Polkinghorne A, Timms P. 2013. Chlamydia pneumoniae: modern insights into an ancient pathogen. Trends Microbiol. 21(3):120–8.

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Mycoplasma pneumoniae

1 – 10 copies of bacterial DNA in the amplified sample

Specimens

Sputum, bronchoalveolar lavage, naso-pharyngeal swab, nasopharyngeal aspirate

Detection technology

Molecular Beacons

Performance study and results

Within the frame of testing the functional characteristics of EliGene® Mycoplasma pneumoniae UNI kit, totally 50 clinical specimens (swabs) were analyzed. The EliGene® Mycoplasma pneumoniae UNI Kit diagnosed correctly 6 specimens as *Mycoplasma pneumoniae* positive. Totally 44 specimens were rightly determined by EliGene® Mycoplasma pneumoniae UNI as *Mycoplasma pneumoniae* negative.

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Mycoplasma hominis, *Mycoplasma genitalium*
1 – 10 copies of bacterial DNA in the amplified sample.

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Molecular Beacons

Performance study and results

Within the frame of testing the functional characteristics of EliGene® Mycoplasma hom/gen UNI kit, totally 60 clinical specimens were analyzed. The EliGene® Mycoplasma hom/gen UNI kit diagnosed correctly as *Mycoplasma hominis* and *Mycoplasma genitalium* positive 7 specimens. Totally 53 specimens were rightly determined by EliGene® Mycoplasma hom/gen UNI kit as *Mycoplasma hominis* and *Mycoplasma genitalium* negative.

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Specificita: 100 %

References / Zdroje

Edberg A, Jurstrand M, Johansson E, Wikander E, Höög A, Ahlqvist T, Falk L, Skov J, Fredlund H. 2008. A comparative study of three different PCR assays for detection of *Mycoplasma genitalium* in urogenital specimens from men and women. Journal of Medical Microbiology, 57, 304–309.

Baczynska A, Svenstrup HF, Fedder J, Birkeland S, Christiansen G. 2004. Development of real-time PCR for detection of *Mycoplasma hominis*. BMC Microbiology, 4:35.

Simpson T, Oh MK. 2004. Urethritis and cervicitis in adolescents. Adolesc Med Clin. 15(2):253-71.

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Chlamydia trachomatis, including "Sweden variant"

1 – 10 copies of bacterial DNA in the amplified sample

Specimens

Urine, swabs (cervical, urethral, conjunctival, rectal), sperm

Detection technology

Molecular Beacons

Performance study and results

110 clinical samples were analysed using the Elisabeth Pharmacon EliGene® Chlamydia trachomatis UNI test kit and the results were verified using the Cobas Amplicor C. trachomatis test. Of the 110 samples, 51 were correctly identified as positive for *Chlamydia trachomatis* and 59 were negative.

Sensitivity: 100%

Specificity: 100%

References / Zdroje

Jurstrand M, Christerson L, Klint M, Fredlund H, Unemo M, Herrmann B. 2010. Characterisation of Chlamydia trachomatis by ompA sequencing and multilocus sequence typing in a Swedish county before and after identification of the new variant. Sex Transm Infect. 86(1):56-60.

Kalwij S, Macintosh M, Baraitser P. 2010. Screening and treatment of Chlamydia trachomatis infections. BMJ. 21:340:c1915.

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Neisseria gonorrhoeae

1 – 10 copies of bacterial DNA in the amplified sample.

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Specimens

Urine, swabs (cervical, urethral, conjunctival, rectal), sperm

Detection technology

Molecular Beacons

Performance study and results

88 Clinical samples were analysed and the results verified using the Cobas Amplicor CT/NG test. The EliGene[®] Neisseria UNI kit correctly identified 35 samples as positive and 53 samples as negative for *Neisseria gonorrhoeae*.

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Ureaplasma urealyticum, *Ureaplasma parvum*
1 – 10 copies of bacterial DNA in the amplified sample.

Specimens

Urine, swabs (cervical, urethral, conjunctival, rectal), sperm

Detection technology

Molecular Beacons

Performance study and results

Within the frame of testing the functional characteristics of EliGene® Ureaplasma UNI kit, totally 68 clinical specimens were analyzed. The EliGene® Ureaplasma UNI Kit diagnosed correctly as *Ureaplasma urealyticum* or *Ureaplasma parvum* positive 20 specimens. Totally 48 specimens were rightly determined by EliGene® Ureaplasma UNI kit as *Ureaplasma urealyticum* or *Ureaplasma parvum* negative.

Sensitivity: 100%

Specificity: 100%

References / Zdroje

Allam AB, Alvarez S, Brown MB, Reyes L. 2011. *Ureaplasma parvum* infection alters filamin A dynamics in host cells. BMC Infect Dis. 11:101.

Pinna GS, Skevaki CL, Kafetzis DA. 2006. The significance of *Ureaplasma urealyticum* as a pathogenic agent in the paediatric population. Curr Opin Infect Dis. 19(3):283-9.

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Specimens

Urine, sputum, bronchoalveolar lavage, exudates, tissue, paraffin embedded tissue

Detection technology

Molecular Beacons

Performance study and results

Within the frame of performance study were tested by EliGene[®] MTB UNI Kit totally 533 clinical samples. EliGene[®] MTB UNI kit detected correctly all 52 positive samples. Totally 481 samples were correctly diagnosed by EliGene[®] MTB UNI Kit as negative.

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
Figure: Sensitivity test of EliGene[®] MTB UNI kit on LightCycler[®] 2.0 instrument shows the detection of 2 DNA copies in reaction mix.

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Borrelia burgdorferi sensu lato

1 – 10 copies of bacterial DNA in the amplified sample

Specimens

Blood, urine, CSF, synovial fluid, tissue

Detection technology

Molecular Beacons

Performance study and results

Within the frame of performance study were tested by EliGene® Borrelia UNI kit 486 clinical samples. EliGene® Borrelia UNI kit detected correctly all 9 positive samples. Totally 477 samples were correctly diagnosed by EliGene® Borrelia UNI kit as *Borrelia* negative.


Sensitivity: 100%

Specificity: 100%

References / Zdroje

Hytönen J, Hartiala P, Oksi J, Viljanen MK. 2008. Borreliosis: recent research, diagnosis, and management. Scand J Rheumatol. 37(3):161-172.

Priem S, Rittig MG, Kamradt T, Burmester GR, Krause A. 1997. An optimized PCR leads to rapid and highly sensitive detection of *Borrelia burgdorferi* in patients with Lyme borreliosis. J Clin Microbiol. 35(3):685–690.

Figures: Detection of pathogenic *Borrelia* sp. in DNA samples (isolated from ticks) by EliGene® Borrelia UNI kit on ABI®7500FAST instrument.

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
Molecular Beacons

Performance study and results

Within the frame of performance study were tested by EliGene[®] Anaplasma phagocytophilum UNI kit 527 samples of ticks. EliGene[®] Anaplasma phagocytophilum UNI kit detected correctly all positive 26 ticks samples. Totally 501 ticks samples were right determined by EliGene[®] Anaplasma phagocytophilum UNI kit as *Anaplasma phagocytophilum* negative.

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Specificita: 100 %

References / Zdroje

Gary P. Wormser et al. 2006. The Clinical Assessment, Treatment, and Prevention of Lyme Disease, Human Granulocytic Anaplasmosis, and Babesiosis: Clinical Practice Guidelines by the Infectious Diseases Society of America. Clin Infect Dis. 43(9):1089-1134.

Johan S. Bakken, MD, Stephen Dumler. 2008. Human Granulocytic Anaplasmosis. Clin Infect Dis. 22(3):433-448

EliGene[®] Anaplasma phagocytophilum UNI
Cat. No.: 90061-UNI 50 reactions

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EliGene® HCV LC

Analytical specificity and sensitivity

Hepatitis C virus
50 IU/ml of serum

Specimens

Serum, plasma

Detection technology

Molecular Beacons

Performance study and results

Within the frame of performance study of EliGene® HCV LC kit overall 500 clinical specimens were analyzed. From these specimens, 300 HCV positive specimens and 200 HCV negative specimens were verified by COBAS AMPLICOR HCV MONITOR Test, version 2.0 (v2.0). The EliGene® HCV LC kit diagnosed as HCV positive 299 specimens and 201 as negative.

100% sensitivity can be declared for the concentration of 50 IU/ml of serum isolated by Magneto 1ml Serum/Plasma DNA/RNA isolation kit - extraction volume: 1.0 ml, elution volume: 50 µl. The linear range of the EliGene® HCV LC kit has been determined to cover concentrations from 1×10^5 IU/ml to at least 50 IU/ml.

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Specificita: 100 %

References / Zdroje

Bannister BA, Begg NT, Gillespie SH. 2000. Infectious Disease. Blackwell Science, 2th Ed.


Figure: Comparison of the COBAS AMPLICOR HCV MONITOR Test (v.2.0) against the EliGene® HCV LC Kit

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
Figure: Calculation of the linear range. The equation of the regression line is included in the figure.

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

Figure: Analysis of specimens with different concentration of HCV from 126 IU/ml to 21000 IU/ml with EliGene® HCV LC kit.
A – amplification curves for five HCV standards,
B – amplification curves for HCV positive specimens.

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EliGene® HBV RT

Analytical specificity and sensitivity

Hepatitis B virus
0,1 IU/ml of isolated specimen

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Molecular Beacons

Performance study and results

Within the frame of performance study of EliGene® HBV RT kit overall 500 clinical specimens were analyzed. From these specimens, 300 HBV positive specimens and 200 HBV negative specimens were confirmed by artus® HBV TM PCR Kit. The EliGene® HBV RT kit diagnosed as HBV positive 299 specimens and 201 as negative.

100% sensitivity can be declared for both kits for the concentration of 5 IU/ml of serum corresponding to concentration of 0,1 IU/ml of isolated specimen by Magneto 1ml Serum/Plasma DNA/RNA isolation kit - extraction volume: 1,0 ml, elution volume: 50 µl. The linear range of the EliGene® HBV RT kit has been determined to cover concentrations from 1×10^5 IU/ml to at least 0,1 IU/ml.

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References / Zdroje

Bannister BA, Begg NT, Gillespie SH. 2000. Infectious Disease. Blackwell Science, 2th Ed.

Ho SKN, Yam W-C, Leung ETK, Wong L-P, Leung JKH, Lai K-N, Chan TM. 2003. Rapid quantification of hepatitis B virus DNA by real-time PCR using fluorescent hybridization probes. J Med Microbiol. 52(5):397-402.



Figure: Calculation of the linear range.

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LightCycler® 2.0 only with EliBlocker

EliGene® Influenza A/B/pandemic LC

Analytical specificity and sensitivity

Influenza virus A, B, A/H1N1

1 – 10 molecules of viral RNA in the amplified sample

Specimens

Serum, plasma, sputum, bronchoalveolar lavage, nasopharyngeal swab / aspirate

Detection technology

Molecular Beacons

Performance study and results

Within the frame of performance study of EliGene® Influenza A/B/Pandemic LC kit overall 100 clinical specimens were analyzed. From these specimens 61 blind specimens and 39 positive specimens were verified by Cepheid Xpert® Flu Kit and INFLUENZA A/B Q-PCR Alert kit. From these 100 samples, 9 samples were Influenza A virus RNA positive, 6 samples Influenza B virus RNA positive and 24 Influenza A/H1N1 virus RNA positive. The EliGene® Influenza A/B/Pandemic LC kit diagnosed as Influenza A, B, A/H1N1 virus RNA positive all 39 and all 61 negative specimens. There were no discrepant results.

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Specificita: 100 %

References / Zdroje

Rose N, Hervé S, Eveno E, Barbier N, Eono F, Dorenlor V, Andraud M, Camsusou C, Madec F, Simon G. 2013. Dynamics of influenza A virus infections in permanently infected pig farms: evidence of recurrent infections, circulation of several swine influenza viruses and reassortment events. Vet Res. 44(1):72.

Labella AM, Merel SE. 2013. Influenza. Med Clin North Am. 97(4):621-45.

Cha RM, Smith D, Shepherd E, Davis CT, Donis R, Nguyen T, Nguyen HD, Do HT, Inui K, Suarez DL, Swayne DE, Pantin-Jackwood M. 2013. Suboptimal protection against highly pathogenic avian influenza (H5N1) viruses from Vietnam in ducks vaccinated with commercial poultry vaccines. Vaccine. S0264-410X(13)01137-7.

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EliGene® HSV1/HSV2 UNI

Analytical specificity and sensitivity

Herpes virus 1 and 2

1 – 10 molecules of viral RNA in the amplified sample

Specimens

CSF, biopsy, swabs (skin, mucosa), blood

Detection technology

TaqMan

Performance study and results

Within the frame of performance study of EliGene® HSV1/HSV2 UNI kit overall 636 HSV1 clinical specimens and 211 HSV2 clinical specimens were analyzed. From 636 HSV1 clinical specimens were 201 positive. EliGene® HSV1/HSV2 UNI kit analyzed 197 samples as positive. Totally 435 clinical specimens were properly by EliGene® HSV1/HSV2 UNI kit analyzed as negative. From 211 HSV2 clinical specimens were 24 positive. EliGene® HSV2 UNI kit analyzed all 24 samples as positive. Totally 187 clinical specimens were properly by EliGene® HSV2 UNI kit analyzed as negative.

HSV1 Sensitivity: 98.04%

HSV1 Specificity: 100%

HSV2 Sensitivity: 100%

HSV2 Specificity: 100%

References / Zdroje

Bannister BA, Begg NT, Gillespie SH. 2000. Infectious Disease. Blackwell Science, 2th Ed.

Hlinomazová Z, Loukotová V, Horáčková M, Šerý O. 2010. The treatment of HSV1 ocular infections using quantitative real-time PCR results. Acta Ophthalmol. 10:1755-3768.

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EliGene® VZV UNI

Analytical specificity and sensitivity

Varicella zoster virus

1 – 10 molecules of viral DNA in the amplified sample

Specimens

CSF, biopsy, swabs (skin, mucosa), blood

Detection technology

TaqMan

Performance study and results

Within the frame of performance study of EliGene VZV UNI kit overall 481 clinical specimens were analyzed. From these 481 clinical specimens 51 were positive. EliGene® VZV UNI kit analyzed all 51 samples as positive. Totally 430 clinical specimens were properly by EliGene® VZV UNI kit analyzed as negative.

Sensitivity: 100%

Specificity: 100%

References / Zdroje

Bannister BA, Begg NT, Gillespie SH. 2000. Infectious Disease. Blackwell Science, 2th Ed.

Dwyer DE, Cunningham AL. 2002. Herpes simplex and varicella-zoster virus infections. Med J Aust. 177(5):267-73.

Hlinomazová Z, Loukotová V, Horáčková M, Šerý O. 2010. The treatment of HSV1 ocular infections using quantitative real-time PCR results. Acta Ophthalmol. 10: 1755-3768.

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EliGene® Adenovirus RT

Analytical specificity and sensitivity

all 51 types of Adenoviruses

1 – 10 molecules of viral DNA in the amplified sample

Specimens

Urine, CSF, feces, swabs, serum, plasma

Detection technology

TaqMan

Performance study and results

Within the frame of performance study of EliGene® Adenovirus UNI kit overall 350 clinical specimens were analyzed. From these 350 clinical specimens were 20 positive. EliGene® Adenovirus UNI kit analyzed 20 samples as positive. Totally 330 clinical specimens were properly by EliGene® Adenovirus UNI kit analyzed as negative.

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EliGene® Norovirus LC

Analytical specificity and sensitivity

Norovirus genogroup 1 and 2

1 – 10 molecules of viral RNA in the amplified sample

Specimens

Serum, plasma, feces, water, CSF

Detection technology

Molecular Beacons

Performance study and results

Within the frame of performance study of EliGene® Norovirus UNI kit overall 100 clinical specimens were analyzed. From these specimens 100 blind specimens were verified by reference method (La Rosa et al. 2010). From these 100 samples 56 samples were Norovirus RNA positive. The EliGene® Norovirus UNI kit diagnosed as Norovirus RNA positive all 56 specimens. There were no discrepant results. Totally 44 specimens were right determined by the EliGene® Norovirus UNI kit as Norovirus RNA negative.

Sensitivity: 100%

Specificity: 100%

References / Zdroje

Bannister BA, Begg NT, Gillespie SH. 2000. Infectious Disease. Blackwell Science, 2th Ed.

Oberste MS, Maher K, Kilpatrick DR, Pallansch MA. 1999. Molecular evolution of the human Noroviruses: correlation of serotype with VP1 sequence and application to picornavirus classification". J. Virol. 73 (3):1941–8

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EliGene[®] Enterovirus LC

Analytical specificity and sensitivity

Enteroviruses (Coxsackie A, Coxsackie B and Echo virus)

1 – 10 molecules of viral RNA in the amplified sample

Analytická specificita a senzitivita

Enteroviruses (Coxsackie A, Coxsackie B and Echo virus)

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Within the frame of performance study of EliGene[®] Enterovirus LC kit overall 100 clinical specimens were analyzed. From these specimens 100 blind specimens were verified by reference method (La Rosa et al. 2010). From these 100 samples 56 samples were Enterovirus RNA positive. The EliGene[®] Enterovirus LC kit diagnosed as Enterovirus RNA positive all 56 specimens. There were no discrepant results. Totally 44 specimens were right determined by the EliGene[®] Enterovirus LC kit as Enterovirus RNA negative.

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Specifita: 100 %

References / Zdroje

Bannister BA, Begg NT, Gillespie SH. 2000. Infectious Disease. Blackwell Science, 2th Ed.

Oberste MS, Maher K, Kilpatrick DR, Pallansch MA. 1999. Molecular evolution of the human enteroviruses: correlation of serotype with VP1 sequence and application to picornavirus classification". J. Virol. 73(3):1941–8.

Santti, Juhana; Heli Harvala, Leena Kinnunen, and Timo Hyypia. 2000. "Molecular epidemiology and evolution of coxsackievirus A9. Journal of General Virology 81(5):1361–1372.


Figure: Analysis of Positive sample by EliGene[®] Enterovirus LC kit on LightCycler[®] 2.0 instrument showing results in 530 channel (Enterovirus detection)

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EliGene® TBEV LC

Analytical specificity and sensitivity

Genome RNA of tick-borne encephalitis virus
50 molecules of viral RNA in the amplified sample

Specimens

Serum, plasma, CSF, ticks,

Detection technology

TaqMan

Pathogen description:

Tick - borne encephalitis (TBE) is a viral infectious disease involving the central nervous system. The disease most often manifests as meningitis, encephalitis, or meningoencephalitis. Although TBE is most commonly recognized as a neurological disorder, mild fever can also occur.

The tick-borne encephalitis virus is known to infect a range of hosts including ruminants, birds, rodents, carnivores, horses and humans. The disease can also be zoonotic.

References / Zdroje

Mansfield, K.L., Johnson, N., Phipps, L.P., Stephenson, J.R., Fooks, A.R., and Solomon, T. (2009). Tick-borne encephalitis virus – a review of an emerging zoonosis. *J Gen Virol* 90, 1781–1794

Lindquist, L., Vapalahti, O. (2008). Tick-borne encephalitis. *Lancet* 371(9627), 1861-71

Suess, J. (2008) Tick-borne encephalitis in Europe and beyond -- the epidemiological situation as of 2007. *Euro Surveillance* 13(26)

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EliGene® Coeliac RT

EliGene® Coeliac RT28

EliGene® Coeliac 3.0 RT is intended for the screening alleles:

HLA-DQ2.5 (DQA1*05/DQB1*02),
HLA-DQ2.2 (DQA1*02/DQB1*02),
HLA-DQ8 (DQA1*03/DQB1*03:02)
from isolated DNA

EliGene® Coeliac RT is intended for the screening alleles:

HLA-DQ2.5 (DQA1*05/DQB1*02),
HLA-DQ8 (DQA1*03/DQB1*03:02),
HLA-DR4 (DRBR1*04)

from isolated DNA

EliGene® Coeliac RT28 is intended for the screening alleles:

HLA-DQ2.5 (DQA1*05/DQB1*02)
HLA-DQ8 (DQA1*03/DQB1*03:02)
from isolated DNA

EliGene® Coeliac 3.0 RT určen pro skrining alel:

HLA-DQ2.5 (DQA1*05/DQB1*02),
HLA-DQ2.2 (DQA1*02/DQB1*02),
HLA-DQ8 (DQA1*03/DQB1*03:02)
z izolované DNA

EliGene® Coeliac RT určen pro skrining alel:

HLA-DQ2.5 (DQA1*05/DQB1*02),
HLA-DQ8 (DQA1*03/DQB1*03:02),
HLA-DR4 (DRBR1*04)

z izolované DNA

EliGene® Coeliac RT28 určen pro skrining alel:

HLA-DQ2.5 (DQA1*05/DQB1*02)
HLA-DQ8 (DQA1*03/DQB1*03:02)
z izolované DNA

Specimens

Buccal swabs, blood

Detection technology

TaqMan LNA

Performance study and results

EliGene® Coeliac RT kit specificity was tested on 50 samples of human DNA with genotypes determined by DNA sequencing. Totally 50 samples were right determined by EliGene® Coeliac RT.

The clinical specificity of EliGene® Coeliac RT kit is 100%.

Analytical sensitivity: from 1 ng DNA molecules
Specificity: 100 %

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Analytická citlivost: od 1 ng DNA molekul
Specifita: 100 %

References / Zdroje

Megiorni F. and Pizzuti A. 2012. HLA-DQA1 and HLA-DQB1 in Celiac disease predisposition: practical implications of the HLA molecular typing. JBS 19: 88-93.

Krini M. Et al. 2012. HLA class II high-resolution genotyping in Greek children with celiac disease and impact on disease susceptibility. Pediatr Res. 72(6):625-30.

Troncone R, Jabri B. 2011. Coeliac disease and gluten sensitivity. J Intern Med. 269(6):582-90.

Trynka G, Wijmenga C, van Heel DA. 2010. A genetic perspective on coeliac disease. Trends Mol Med. 6(11):537-50.

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EliGene® Spondylitis HLA-B27 RT

EliGene® Spondylitis HLA-B27 C6

Analytical specificity and sensitivity
allele HLA-B27

Analytická specificita a senzitivita
alela HLA-B27

Specimens

Buccal swabs, blood

Vzorky

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Detection technology

TaqMan

Technologie detekce

TaqMan

Performance study and results

Within the frame of performance study of EliGene® HLA-B27 RT kit overall 226 specimens were analyzed. From these specimens 226 blind specimens were verified by reference method (Roelandse-Koop EA et al., 2011). From these 226 samples 36 samples were HLA-B27 allele positive. The EliGene® HLA-B27 RT kit diagnosed as HLA-B27 allele positive positive all 36 specimens and 190 specimens were right determined as HLA-B27 allele negative.

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EliGene® Spondylitis HLA-B27 RT / C6

Cat. No.: 90060-RT	50 reactions
Cat. No.: 90060-C6	96 reactions

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
C6 version: Cobas® 6800/8800 Systems

ONCOGENETICS

EliGene® Colorectum NGS	33
EliGene® Lung NGS	34
EliGene® GIST NGS	35
EliGene® Adaptor-IL NGS (accessories)	36

ONKOGENETIKA

EliGene® Colorectum NGS

Analytical specificity and sensitivity

Kit is intended for multiplex PCR of target regions of **KRAS** (exons 2,3 and 4), **NRAS** (exons 2,3 and 4) and **BRAF** genes (exon 15) allowing detection of polymorphisms included in the amplified regions.

Analytical sensitivity is 20 ng of DNA in reaction mix.

Specimens

Blood, FFPE material

Detection technology

Multiplex PCR amplification of all required target regions of the gene(s) of interest with consequent massively parallel sequencing (MPS) using NGS instrument Illumina MiSeq.

Performance study and results

Overall 12 isolates of FFPE tissue were analysed by **EliGene® Colorectum NGS** kit together with **EliGene® Adaptor-IL NGS** kit for the presence of clinically significant mutations in target regions of KRAS (exons 2,3 and 4), NRAS (exons 2,3 and 4) and BRAF genes (exon 15) within the frame of performance study. These were compared with previous results obtained by reference methods KRAS StripAssay (ViennaLab), Cobas KRAS Mutation Test (Roche), NGS-Somatic1 (Multiplicom), LightMix KRAS NRAS (TIBMolBiol), Cobas EGFR Mutation Test (Roche) and BRAF StripAssay (ViennaLab). The same mutations were found as in previously analysed samples. The number of obtained reads was verified to be higher than 100 read pairs with the mean number of 3000 read pairs/detected mutation. The reproducibility of testing was proved in two independent sequencing runs with different lots.

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EliGene[®] Lung NGS

Analytical specificity and sensitivity

Kit is intended for multiplex PCR of target regions of EGFR gene (exons 18, 19, 20, 21) allowing detection of polymorphisms included in the amplified regions.

Analytical sensitivity is 20 ng of DNA in reaction mix.

Specimens

Blood, FFPE material

Detection technology

Multiplex PCR amplification of all required target regions of the gene(s) of interest with consequent massively parallel sequencing (MPS) using NGS instrument Illumina MiSeq.

Performance study and results

Overall 12 isolates of FFPE tissue were analysed by EliGene[®] Lung NGS kit together with EliGene[®] Adaptor-IL NGS kit for the presence of clinically significant mutations in target regions of EGFR gene (exons 18, 19, 20, 21) within the frame of performance study. These were compared with previous results obtained by reference methods Cobas EGFR Mutation Test (Roche) and Sanger sequencing. The same mutations were found as in previously analysed samples. The number of obtained reads was verified to be higher than 100 read pairs with the mean number of 3000 read pairs/detected mutation. The reproducibility of testing was proved in two independent sequencing runs with different lots.

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EliGene® GIST NGS

Analytical specificity and sensitivity

Kit is intended for multiplex PCR of target regions of genes **c-KIT** (exons 9, 11, 13, 14, 15, 16, 17) and **PDGFR α** (exons 8, 10, 12, 14, 18) allowing detection of polymorphisms included in the amplified regions.

Analytical sensitivity is 20 ng of DNA in reaction mix.

Specimens

Blood, FFPE material

Detection technology

Multiplex PCR amplification of all required target regions of the gene(s) of interest with consequent massively parallel sequencing (MPS) using NGS instrument Illumina MiSeq.

Performance study and results

Overall 12 isolates of FFPE tissue were analysed by EliGene® GIST NGS kit together with EliGene® Adaptor-IL NGS kit for the presence of clinically significant mutations in target regions of c-Kit (exons 9, 11, 13, 14, 15, 16, 17) and PDGFR α genes (exons 8, 10, 12, 14, 18) within the frame of performance study. These were compared with previous results obtained by reference methods NGS-Somatic1 (Multiplicom), Cobas EGFR Mutation Test (Roche) and Sanger sequencing. The same mutations were found as in previously analysed samples. The number of obtained reads was verified to be higher than 100 read pairs with the mean number of 3000 read pairs/detected mutation. The reproducibility of testing was proved in two independent sequencing runs with different lots.

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EliGene[®] Adaptor-IL NGS

Intended use

Kit is intended for the incorporation of molecular identifiers (MIDs) and specific sequencing adaptors to all amplicons generated using the multiplex PCR based on EliGene[®] NGS kits, enabling the identification of polymorphisms in the target area.

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Specimens

Suitable specimens are the amplicons from multiplex PCR using EliGene[®] Lung NGS, EliGene[®] Colorectum NGS and EliGene[®] GIST NGS.

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Illumina MiSeq, MiniSeq

2. Product guarantee policy

Product quality is our highest priority at **ELISABETH PHARMACON**. Clinical studies are available for all CE certified kits. As manufacturing, testing, packaging and delivery are carried out under the **EN ISO 13485:2016** quality system, we guarantee batch-to-batch consistency. All **EliGene** products are supplied with a Certificate of Analysis, summarizing the conditions and results of quality assurance tests.

3. Shipping

EliGene diagnostic kits must be shipped frozen on dry ice, for short journeys, kits can be shipped packed on blue ice. Freeze-thaw cycles lower the sensitivity of diagnostic kits. Diagnostic kits are packed in polystyrene boxes to protect them from thawing during shipment and laboratory work.

4. Advantages of EliGene® detection kits

- The usage of the Internal Control allows to monitor the extraction procedure and to check possible PCR inhibition
- Kits are supplied with positive control allowing quality control of the analysis
- Simultaneous amplification (multiplex) of DNA/RNA from pathogen and internal control in one PCR tube
- The usage of Hot-Start technology minimize the risk of non-specific reaction and provide ultimate sensitivity 1-10 copies of nucleic acid/reaction
- The usage of Ready to use Mastermixes containing all reaction components for easier reaction setup and handling
- The usage of uracil in the Mastermixes of UNI kits enables UNG treatment to prevent the re-amplification of carryover-PCR products by removing any uracil incorporated into single- or double-stranded DNA
- To ensure the highest quality of our kits, the raw material from the leading providers is used
- Quantitative Standards in the most of the kits enables to calculate precisely the NA load from 10² to 10⁶ copies/reaction
- High stability of used reagents
- The usage of EliGene® UNI kits enables complete detection in 60 minutes
- Kits are CE certified IVD
- Kits are validated on LightCycler® 2.0, Nano and 480, ABI® 7300, and 7500FAST, RotorGene-Q(RotorGene 6000) and MyGo® cycler.
- Compatibility of the kits with broad range of instruments - RotorGene®, SmartCycler®, MX3000P® and 3005P® QPCR systems or iQ5®, CFX96® and other cyclers.

LightCycler is registered trademark of Roche Molecular Systems Inc., "ABI" is a registered trademark of Applied Biosystems. "iQ5", "MyiQ" and "CFX96 & CFX384" are registered trademarks of Bio-Rad Laboratories. "SmartCycler" is registered trademark of Cepheid. "Rotor-Gene" is a trademark of Corbett Research. "MX 3000P" and "MX3000P" are registered trademarks of Stratagene, Inc. MyGo® is registered trademark of IT-IS Life Science Ltd.

5. Sample processing

The amplification success depends mostly on the quality of the purified NA. The highest sensitivity in PCR detection can only be achieved when you have first-rate sample preparation. The quality of sample preparation depends not only on isolation procedure but also on timing of collection (clinical relevancy), transport and storage conditions. The presence of conservation agents in transport media could lead to PCR inhibition as well. Taking together all these facts, the isolation procedure must be rigorously tested to avoid false negative results.

7. Product lines of EliGene kits

EliGene® kits for Real-Time PCR


EliGene® kits for real-time PCR analysis are divided according to detection chemistry to basic lines RT, UNI and LC virus, respectively.

In **RT line**, the TaqMan chemistry using hydrolyzing fluorescence probes is used in the combination with hot-start polymerase and ROX dye as a passive reference (see Figure 1A). The total reaction volume is 25ul when 5ul of isolated DNA is added directly to ready to use mastermix containing all necessary components for reaction. In RT kits, the fluorescence dyes FAM and HEX with excitation wavelengths of 495, 520 nm and emission wavelengths of 520, 550 nm together with non-fluorescence quencher BHQ1 are used. These kits are basically designed for ABI instruments (line 7000 and SteOne) and they are optimized on LightCycler® 480 instrument, however can be used with the other instruments that meet the spectral characteristics of used fluorescence dyes. The typical time of the run for RT kits is about 105 minutes.

In **UNI line**, the Molecular Beacons using hybridization fluorescence probes or TaqMan chemistry using hydrolyzing fluorescence probes are used in combination with hot-start polymerase (see Figure 1). The total reaction volume is 20 ul when 5 ul of isolated DNA is added directly to ready-to use mastermix containing all necessary components for reaction. The usage of uracil in the reaction mix enables UNG treatment to prevent the re-amplification of carryover-PCR products by removing any uracil incorporated into single- or double-stranded DNA. In UNI kits the fluorescence dyes FAM and HEX with excitation wavelengths of 495, 520 and emission wavelengths of 520, 550 together with non-fluorescence quencher BHQ1 are used. These kits are designed for almost all real-time instruments on the market. The typical time of run for UNI kits is about 60 minutes.

In **LC virus line**, the Molecular Beacons using hybridization fluorescence probes or TaqMan chemistry using hydrolyzing fluorescence probes are used in combination with one-step master mix for reverse transcription with consequent qPCR analysis (see Figure 1). The total reaction volume is 20/30 ul when 5/10 ul of isolated RNA is added to mastermix created by mixing of polymerase with reaction mix containing all necessary components for the reaction. The fluorescence dyes FAM, HEX and Cy5 with excitation wavelengths of 495, 520 and 649 and emission wavelengths of 520, 550 and 666 together with non-fluorescence quenchers BHQ1 and BHQ2 are used. The LC kits. These kits are designed for the majority of real-time instruments on the market. The typical time of run for LC kits is about 70 minutes.

Figure 1. Basic principle of TaqMan hydrolyzation probes (A) and Molecular Beacons hybridization probes (B) action.

The TaqMan probe hybridizes between the primers in annealing step and during elongation step probe is cleaved by exonuclease activity of Taq Polymerase. The Molecular Beacon probe forms beacon-like structure that opens to specific sequence during the hybridization in annealing step (B).

"TaqMan" is registered trademark of Roche Molecular Systems Inc., "Molecular Beacon" is registered trademark of Public Health Research Institute Properties, Inc. (PHRI)

EliGene® kits for Next Generation Sequencing

EliGene® NGS kits are used for one or more multiplex PCR reactions resulting in amplification of target regions of one sample (Figure 2A). The resulting amplicons are barcoded with MID sequences, and p7, p5 adaptors required for Illumina MiSeq Sequencing (Figure 2B). Each amplicon library is subsequently purified from small residual DNA fragments shorter than 150bp and the DNA concentration is determined using a Real-Time PCR. Next, these purified and individually tagged amplicon libraries are pooled equimolar, resulting in an amplicon pool or sequencing sample, which is sequenced using a NGS instrument Illumina for massively parallel sequencing (MPS) according to the manufacturer's instructions. The resulting sequence reads are subsequently analysed by comparison with the reference sequence of the targeted gene(s) to identify variant positions. For amplification is used a ready-to-use mastermix containing a hot-start DNA polymerase. The positions of the individual custom sequencing primers supplied by EliGene® NGS kit as well as the Illumina sequencing reagents are indicated and clarified in Figure 2C.


Figure 2. Layout of the process resulting in multiplex PCR of target regions (A), incorporation of MIDs and p7 and p5 adaptors (B) and structure of resulting product after the sequencing on Illumina MiSeq instrument (C)

Product portfolio of ELISABETH PHARMACON

Portfolio of **EliZyme™** enzymes:

End-Point PCR

- HotStart Robust
- HotStart FAST
- FAST Taq
- HIFI
- OneS

Real-Time PCR

- Green MIX
- Probe MIX
- OneS Green
- OneS Probe
- HRM
- Genotyping

NEXT GENERATION SEQUENCING

- Library Quantification Kit
- Hyper Preparation Kit
- HIFI HotStart MIX

Portfolio polymeráz **EliZyme™**

End-Point PCR


- HotStart Robust
- HotStart FAST
- FAST Taq
- HIFI
- OneS

Real-Time PCR

- Green MIX
- Probe MIX
- OneS Green
- OneS Probe
- HRM
- Genotyping

NEXT GENERATION SEQUENCING

- Library Quantification Kit
- Hyper Preparation Kit
- HIFI HotStart MIX

EliGene® Isolation kits are available in the version for research only or certified for clinical use. In our portfolio there are kits for manual DNA and RNA isolation on column from different materials.

PORTFOLIO OF KITS

- EliGene® MTB Isolation Kit (CE)
- EliGene® Urine Isolation Kit (CE)
- EliGene® Soil DNA Isolation Kit
- EliGene® Blood DNA Isolation Kit
- EliGene® Plant DNA Isolation Kit
- EliGene® Tissue DNA Isolation Kit
- EliGene® FFPE Tissue DNA Isolation Kit
- EliGene® Microbial DNA Isolation Kit
- EliGene® Viral DNA/RNA Isolation Kit
- EliGene® BodyFluid DNA Isolation Kit

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PORTFOLIO SOUPRAV

- EliGene® MTB Isolation Kit (CE)
- EliGene® Urine Isolation Kit (CE)
- EliGene® Soil DNA Isolation Kit
- EliGene® Blood DNA Isolation Kit
- EliGene® Plant DNA Isolation Kit
- EliGene® Tissue DNA Isolation Kit
- EliGene® FFPE Tissue DNA Isolation Kit
- EliGene® Microbial DNA Isolation Kit
- EliGene® Viral DNA/RNA Isolation Kit
- EliGene® BodyFluid DNA Isolation Kit

Certificates ELISABETH PHARMACON


ELISABETH PHARMACON® has certificated manufacturing and R&D department by ISO 9001:2015 and ISO 13485:2016.

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Elisabeth Pharmacon CZ
Elisabeth Pharmacon UK
Elisabeth Pharmacon SK
members of Elisabeth Pharmacon Group

www.elisabethgroup.com | www.elisabeth.cz | www.eligene.com | www.elizyme.com

Manufacturer:

ELISABETH PHARMACON Ltd.
Rokycanova 4437, 615 00 Brno-Zidenice
phone: +420 542 213 851
e-mail: info@elisabeth.cz
www.elisabeth.cz

Distributor:

EUROPEAN UNION
European Regional Development Fund
Operational Programme Enterprise
and Innovations for Competitiveness