PRODUCT CATALOG POLYMERASES

End-Point PCR
Real-Time PCR
NGS

Introduction

Dear customers,

ELISABETH PHARMACON is innovative biotechnological company developing modern products for clinical and scientific laboratories. In our portfolio you can find: kits for DNA diagnostics, human genetics, DNA and RNA isolation and plenty of special polymerases – these are presented in this catalogue.

Thanks to positive feedback and affirmative response from satisfied customers, we dare to say that EliZyme polymerases are one of the best polymerases in the world market.

There are end-point PCR kits for difficult DNA templates or long templates up to 35 kb (EliZyme HS Robust), ultra-fast polymerase thermocyclers (EliZyme HS FAST or EliZyme FAST Taq) presented in this catalogue. These enzymes reach speeds up to 1s/kb. Most polymerases are designed as Hot-Start enzymes.

In our portfolio you can find also customer's favourite polymerase with proofreading activity – EliZyme HIFI. Of course, MIX and MIX RED versions are also available. They are making the work with EliZyme polymerases more comfortable.

Polymerases for Real-Time PCR applications are manufactured as mixes. Passive reference dye ROX is included separately making the polymerase compatible on all types of instruments. In the portfolio you can find polymerases with non-PCR inhibiting intercalating dye, for work with probes and enzyme assay for HRM analysis. We also included the one-step reverse transcription kits called OneS.

Whole range of new generation enzymes EliZyme mentioned in this catalogue will surely cover all your needs for various PCR applications.

Content

End-Point PCR Kits	3
Real-Time PCR Kits	11
Next Generation Sequencing Kits	19
Accessories	23

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Obsah

EUROPEAN UNION
European Regional Development Fund
Operational Programme Enterprise
and Innovations for Competitiveness

END-POINT PCR

EliZyme FAST Taq

EliZyme FAST Taq is a robust enzyme suited for everyday PCR applications including amplification of DNA for Sanger sequencing and other genotyping applications. EliZyme FAST Taq is characterised by enhanced PCR speed, yield and specificity. The combination of enzyme and buffer system makes exceptional PCR performance on complex templates such as mammalian genomic DNA.

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5

EliZyme HS FAST

EliZyme HS FAST uses advanced Hot-Start technology for greater sensitivity. Primer-dimer formation and non-specific amplification from low copy number target sequences are avoided by inactivation below 65° C. For higher comfort is EliZyme HS FAST DNA polymerase also available in ready mix version. Ready mix red contains a red dye for tracking during agarose gel electrophoresis. It is suitable for direct loading onto agarose gel.

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6

EliZyme HS Robust

EliZyme HS Robust enables amplification of templates that are GC/AT rich, in low abundance or with PCR inhibitors. Amplification is achievable up to 35kb. For higher comfort is EliZyme HS Robust DNA polymerase also available in a ready mix version. Ready mix red contains a red dye for tracking during agarose gel electrophoresis. It is suitable for direct loading onto agarose gel.

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7

EliZyme OneS Kit

EliZyme OneS Kit includes components for efficient reverse transcription and PCR in a single tube. The PCR is performed by EliZyme HS FAST. EliZyme HS FAST uses advanced Hot-Start technology for greater sensitivity. Primer-dimer formation and non-specific amplification from low copy number target sequences are avoided by inactivation below 65° C.

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8

EliZyme HIFI

EliZyme HIFI is characterised by its 3'-5' exonuclease (proofreading) activity in PCR. The polymerase was derived from Pfu DNA polymerase. The error rate of EliZyme HIFI is 100x lower than of Taq DNA polymerase.

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9

EliZyme FAST Taq

EliZyme FAST Taq is a robust enzyme system suited for routine PCR, amplification of DNA for Sanger sequencing and other genotyping applications.

The enzyme system is characterised by enhanced PCR speed, yield and specificity. The buffer and ready mixes include all components for a successful DNA amplification. EliZyme FAST Taq delivers exceptional PCR performance on complex templates, including GC-rich and AT-rich sequences. EliZyme FAST Taq DNA Polymerase has 5'-3' exonuclease activity but no 3'-5' exonuclease (proofreading) activity. The error rate is approximately 1 error per 2.0×10^5 nucleotides incorporated. PCR products generated with EliZyme FAST Taq are A-tailed and may be cloned into TA vectors.

For higher comfort is EliZyme FAST Taq also available as a 2x ready mix and as EliZyme FAST Taq (no dNTP). EliZyme FAST Taq MIX Red contains a red dye for tracking during agarose gel electrophoresis. It is suitable for direct loading onto agarose gel.

Advantages

- Higher speed (with speed 1s/kb)
- Higher yields under standard and fast PCR conditions
- Increased PCR success rates

Applications

- Routine PCR
- Standard and fast PCR
- Amplification from complex templates including GC-rich and AT-rich sequences
- Sanger sequencing
- TA cloning

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200ng 66ng 22ng 7.3ng 2.4ng 800pg 266pg

Figure

Amplification of a 60% 1.1 kb fragment of GC GAPDH from human genomic DNA. A 3 fold dilution series of template starting from 200ng was used. The final dilution is 266pg. EliZyme FAST Taq is able to amplify lower concentration DNA template compared with competitor "C" and competitor "B".

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EliZyme HS FAST

EliZyme HS FAST uses Hot-Start technology to inactivate the enzyme below 65° C preventing primer-dimer formation and non-specific amplification. DNA polymerase in EliZyme HS FAST is inactivated until the initial activation step at 95° C.

EliZyme HS FAST is suited for difficult PCR templates. The mix is resistant to PCR inhibitors allowing direct PCR from unprocessed samples including bacterial culture, bacterial colonies, blood and urine. EliZyme HS FAST is a robust enzyme system suited for routine PCR, multiplex PCR, amplification of DNA for Sanger sequencing and other genotyping applications.

The enzyme system is characterised by enhanced PCR speed, yield and specificity. The buffer and ready mixes include all components for a successful DNA amplification. EliZyme HS FAST delivers exceptional PCR performance on complex templates including GC-rich and AT-rich sequences. EliZyme HS FAST Taq DNA Polymerase has 5'-3' exonuclease activity but no 3'-5' exonuclease (proofreading) activity. The error rate is approximately 1 error per 2.0×10^5 nucleotides incorporated. PCR products generated with EliZyme HS FAST are A-tailed and may be cloned into TA vectors.

For higher comfort is EliZyme HS FAST also available as a 2x ready mix. EliZyme HS FAST MIX Red contains a red dye for tracking during agarose gel electrophoresis. It is suitable for direct loading onto agarose gel.

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Advantages

- Hot-Start technologie
- Enzyme inactivation below 65° C
- Direct PCR
- Higher speed (with speed 1s/kb)
- Higher yields under standard and fast PCR conditions
- Increased PCR success rates
- Inhibitor tolerant PCR

Applications

- Routine and multiplex PCR
- Direct PCR from bacterial culture, blood and urine
- Amplification from complex templates including GC-rich and AT-rich sequences
- Sanger sequencing
- TA cloning

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
Figure

Amplification of a 1.2 kb fragment of Beta-Actin under standard cycling conditions from human genomic DNA. Hot-Start technology prevents primer-dimer formation. Primer extension redirects DNA polymerase activity from the amplicon of interest leading to reduced sensitivity in the PCR reaction. A 10 fold dilution series of template starting from 100 ng was used. EliZyme HS FAST is able to amplify lower concentration DNA template compared with competitor "A".

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100ng 10ng 1ng 100pg 10pg 1pg 0.1pg 0.01pg

EliZyme HS Robust

EliZyme HS Robust is developed for the amplification of extremely difficult templates. Exceptional performance is delivered by higher processivity and Hot-Start technology. EliZyme HS Robust enables amplification of templates that are GC-rich and AT-rich or with PCR inhibitors. Amplification is achievable up to 35 kb.

The enzyme is inactivated below 65°C preventing primer-dimer formation and non-specific amplification. DNA polymerase in EliZyme HS Robust is inactivated until the initial activation step at 95°C.

The enzyme system is characterised by enhanced PCR speed, yield and specificity. The buffer and ready mixes include all components for a successful DNA amplification. EliZyme HS Robust delivers exceptional PCR performance on complex templates including GC-rich and AT-rich sequences. EliZyme HS Robust Polymerase has error rate approximately 1 error per 5.0×10^5 nucleotides incorporated. PCR products generated with EliZyme HS Robust are A-tailed and may be cloned into TA vectors.

For higher comfort is EliZyme HS Robust also available as a 2x ready mix. EliZyme HS Robust MIX Red contains a red dye for tracking during agarose gel electrophoresis. It is suitable for direct loading onto agarose gel.

Advantages

- Hot-Start technology
- Enzyme inactivation below 65°C
- Increased PCR success rates with amplicons up to 35kb
- Inhibitor tolerant PCR from crude samples
- Higher yields under standard and fast PCR conditions

Applications

- Long Range PCR
- Extremely difficult templates
- Crude sample PCR
- Amplification from complex templates
- Sanger sequencing
- TA cloning

200ng 100ng 50ng 25ng 12.5ng

Figure

Amplification of a 27kb fragment of the p53 gene region from human genomic DNA. A 2 fold dilution series of template starting from 200ng was used. EliZyme HS Robust detects as low as 12.5ng, which is lower than equivalent products from competitors "D" and "B".

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
EliZyme OneS Kit

EliZyme OneS Kit includes components for efficient reverse transcription and PCR in a single tube. MMLV reverse transcriptase (RTase) is thermostable, highly active and has RNase inhibitor. The RTase is not inhibited by ribosomal and transfer RNAs.

The PCR is performed by EliZyme HS FAST Taq DNA Polymerase. The enzyme uses Hot-Start technology to inactivate the enzyme below 65 °C preventing primer-dimer formation and non-specific amplification from target sequences. The enzyme is inactivated until the initial activation step at 95 °C.

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Advantages

- cDNA synthesis and PCR in a single tube
- RNase inhibitor
- Includes EliZyme HS FAST
- Hot-Start technology
- Enzyme inactivation below 65 °C

Applications

- Reverse transcription
- Routine PCR
- Amplification from complex templates including GC-rich and AT-rich sequences
- Low copy number target genes
- Sanger sequencing

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
Figure

Amplification of a 1.2kb fragment of Beta-Actin under standard cycling conditions from human genomic DNA. Hot-Start technology prevents primer-dimer formation. Primer extension redirects DNA polymerase activity from the amplicon of interest leading to reduced sensitivity in the PCR reaction. A 10 fold dilution series of template starting from 100ng was used. EliZyme HS FAST is able to amplify lower concentration DNA template compared with competitor "A".

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100ng 10ng 1ng 100pg 10pg 1pg 0.1pg 0.01pg

EliZyme HIFI

EliZyme HIFI is characterised by its 3'-5' exonuclease (proofreading) activity in PCR. Significantly improved performance is attributed to several point-mutations improving the performance compared with its native form.

EliZyme HIFI is a robust enzyme system suited for routine PCR, DNA amplification for Sanger sequencing and other genotyping applications. The enzyme system is characterised by enhanced PCR speed, yield and specificity. EliZyme HIFI delivers exceptional PCR performance on complex templates including GC-rich and AT-rich sequences. The error rate of EliZyme HIFI is lower than Taq DNA polymerase (1 error per 4.5×10^7 nucleotides incorporated).

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Advantages

- Higher fidelity
- Proofreading activity
- Increased PCR success rates
- Higher yields

Applications

- High fidelity PCR
- Standard and fast PCR
- Amplification from complex templates including GC-rich and AT-rich sequences
- Sanger sequencing

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100ng 50ng 25ng 12.5ng 6.3ng 3.2ng


Figure

Amplification of a 4.8kb fragment of GAPDH from human genomic DNA. A 2 fold dilution series of template starting from 100 ng was used. Cycling conditions: initial denaturation at 90°C for 1 minute, 25 cycles of denaturation at 95°C for 30 seconds, annealing at 60°C for 30 seconds and extension at 72°C for 75 seconds. EliZyme HIFI is able to amplify lower concentration DNA template compared with competitors "B" and "C".

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REAL-TIME PCR

EliZyme Green MIX AddROX

Intercalating dye used in EliZyme Green MIX AddROX does not inhibit PCR, unlike other well-known fluorescent dyes including SYBR® Green. EliZyme Green MIX AddROX combines advanced enzyme, Hot-Start and reaction buffer technology.

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13

EliZyme Probe MIX AddROX

EliZyme Probe MIX AddROX is a universal probe kit designed for use in all probe based Real-Time PCR assays. It is suitable for a multiplex assay or genotyping experiment.

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14

EliZyme OneS Green LowROX

EliZyme OneS Green LowROX Kit is used for quantification of mRNA, total RNA, viral sequences and it specifically detects extremely low copy number targets with high efficiency.

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15

EliZyme OneS Probe Kit

EliZyme OneS Probe Kit is used for quantification of mRNA, total RNA, viral sequences and it specifically detects extremely low copy number targets with high efficiency.

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16

EliZyme HRM MIX

EliZyme HRM MIX uses intercalating dye that does not inhibit PCR. It is essential for accurate SNP genotyping and CpG methylation analysis.

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17

EliZyme Genotyping MIX

EliZyme Genotyping MIX is used in dual-labeled probe based genotyping assays including TaqMan®, molecular beacon® and Scorpion® probe genotyping assays.

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18

EliZyme Green MIX AddROX

Intercalating dye used in EliZyme Green MIX AddROX does not inhibit Real-Time PCR. Primer-dimer formation and non-specific amplification are avoided by Hot-Start technology. EliZyme Green MIX AddROX can be used to quantify any DNA template.

EliZyme Green MIX AddROX is compatible on all Real-Time PCR platforms (under standard and fast cycling conditions).

For higher comfort is mix without ROX. The ROX dye is included separately with the mix and can be added according to your needs.

Advantages

- Non-PCR inhibiting intercalating dye for better signal
- Rapid extension rate for early Ct values
- Increased limit of detection
- Compatible with all Real-Time PCR platforms - standard and fast cycling conditions

Applications

- Absolute quantification
- Relative gene expression analysis
- High-throughput Real-Time PCR from genomic DNA, cDNA and viral sequences
- Low copy number target genes
- Amplification from complex templates

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

Figure

Amplification and melt analysis of PGK1 mouse housekeeping gene. A 3 fold dilution series after extra preincubation at 50 °C was used. Fast cycling conditions: extra preincubation at 50 °C for 10 minutes, initial denaturation at 95 °C for 3 minutes, 45 cycles of denaturation at 95 °C for 3 seconds and annealing/extension at 60 °C for 10 s. Under fast cycling conditions EliZyme Green MIX AddROX shows better sensitivity than competitor "A" and less primer-dimer formation, indicating better Hot-Start technology. Black – competitor "A", Green – EliZyme Green MIX AddROX.

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EliZyme Probe MIX AddROX

EliZyme Probe MIX AddROX is designed for all probe based Real-Time PCR assays including TaqMan®, molecular beacons® and Scorpions® probes. EliZyme Probe MIX AddROX is suitable for a multiplex assay or genotyping experiment. It allows efficient amplification of GC-rich and AT-rich templates. Primer-dimer formation and non-specific amplification are avoided by inhibitor technology.

EliZyme Probe MIX AddROX is compatible on all Real-Time PCR platforms (under standard and fast cycling conditions).

For higher comfort is mix without ROX. The ROX dye is included separately with the mix and can be added according to your needs.

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Advantages

- High efficiency in multiplex reactions
- Rapid extension rate for early Ct values
- Increased limit of detection
- Amplification from complex templates including GC-rich and AT-rich sequences
- Compatible with all Real-Time PCR platforms - standard and fast cycling conditions

Applications

- Absolute quantification
- Relative gene expression analysis
- TaqMan®, Scorpions® and molecular beacon® probes
- Low copy number target genes
- Multiplex or singleplex

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
Figure

The gActin gene was used for singleplex experiments and progress was measured by fluorescence probe. cDNA created from mouse liver RNA was used as template at three concentrations, with three replicates at each concentration. EliZyme Probe MIX AddROX displayed higher fluorescence intensity. Cycling conditions: initial denaturation at 95 °C for 2 minutes, 45 cycles of denaturation at 95 °C for 5 seconds and annealing/extension at 60 °C for 30 s. Black – competitor "A", Green – EliZyme Probe MIX AddROX.

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EliZyme OneS Green LowROX

EliZyme OneS Green LowROX Kit includes components for both cDNA synthesis and Real-Time PCR in a single tube. MMLV reverse transcriptase (RTase) is thermostable, extremely active and has RNase inhibitor. EliZyme OneS Green LowROX kit is used for quantification of mRNA, total RNA and viral sequences. The RTase is not inhibited by ribosomal and transfer RNAs.

Intercalating dye used in EliZyme OneS Green LowROX Kit does not inhibit PCR. Primer-dimer formation and non-specific amplification are avoided by antibody Hot-Start technology.

EliZyme OneS Green LowROX kit is compatible on Real-Time PCR platforms that require low ROX concentration (under fast and ultra fast cycling conditions).

Advantages

- Thermostable reverse transcriptase 45°C to 55°C
- RNase inhibitor
- Non-PCR inhibiting intercalating dye
- Rapid extension rate for early Ct values
- Increased limit of detection

Applications

- Absolute quantification
- Relative gene expression analysis
- Low copy number target genes
- Recommended for template amounts of 10pg – 100ng total RNA or > 0.01pg mRNA per reaction


Figure
Comparison of EliZyme OneS Green LowROX and competitor "B". Amplification of the ACTG1 gene from a dilution series of total RNA extracted from mouse liver. Total RNA concentration varied from 25pg to 250ng. Cycling conditions: reverse transcription at 45°C for 10 minutes, initial denaturation at 95°C for 2 minutes, 45 cycles of denaturation at 95°C for 10 seconds and annealing/extension at 60°C for 10 seconds on Roche LC480. EliZyme OneS Green LowROX had equal performance at high RNA concentrations and superior performance at lower RNA concentrations and lower prevalence of primer-dimer formation. Blue – competitor "B", Green – EliZyme OneS Green LowROX.

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
EliZyme OneS Probe Kit

EliZyme OneS Probe Kit includes components for efficient reverse transcription and Real-Time PCR in a single tube. MMLV reverse transcriptase (RTase) is thermostable, highly active and has RNase inhibitor. EliZyme OneS Probe Kit is used for quantification of total RNA, mRNA, and viral sequences. The RTase is not inhibited by ribosomal and transfer RNAs.

EliZyme OneS Probe Kit is designed for all probe based Real-Time PCR assays including TaqMan[®], molecular beacons[®] and Scorpions[®] probes. Primer-dimer formation and non-specific amplification are avoided by Hot-Start technology.

EliZyme OneS Probe Kit is compatible on Real-Time PCR platforms that do not require addition of ROX as passive reference dye, under standard and fast cycling conditions.

Advantages

- Thermostable reverse transcriptase 45°C to 55°C
- RNase inhibitor
- Rapid extension rate for early Ct values
- Increased limit of detection

Applications

- Absolute quantification
- Relative gene expression analysis
- Low copy number target genes
- TaqMan[®], Scorpions[®] and molecular beacon[®] probes
- Multiplex or singleplex

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

EliZyme HRM MIX

High Resolution Melting (HRM) analysis is a powerful technique used for detection of DNA sequence variations in double-stranded DNA samples.

Intercalating dye used in EliZyme HRM MIX does not inhibit Real-Time PCR. It is essential for accurate SNP genotyping and CpG methylation analysis. The dye has the same affinity for AT-rich and GC-rich regions. Even a difference between A to T can be accurately genotyped.

EliZyme HRM MIX is compatible on particular Real-Time PCR platforms (under standard and fast cycling conditions).

Advantages

- Hot-Start technology for better sensitivity
- Standard and fast cycling conditions

Applications

- Accurate SNP genotyping
- Gene scanning
- Amplification from complex templates including GC-rich and AT-rich sequences

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
Figure

Three samples of human genomic DNA were amplified using EliZyme HRM MIX with primers specific for a fragment of beta-globulin gene. The primers flanked a class IV single nucleotide polymorphism (17A/T). The curves show clear allele calling of a class IV SNP. The black curves are for homozygous A, the purple curves for homozygous T and the green curves for heterozygous A/T.

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EliZyme Genotyping MIX

EliZyme Genotyping MIX is designed for dual-labeled probes including LNA and MGB probes. The mix is used for single nucleotide polymorphism (SNP) genotyping. EliZyme Genotyping MIX is highly efficient in multiplexed reactions. It allows efficient amplification of GC-rich and AT-rich sequences. Primer-dimer formation and non-specific amplification are avoided by inhibitor technology.

EliZyme Genotyping MIX is compatible on all Real-Time PCR platforms (under standard and fast cycling conditions).

EliZyme Genotyping MIX is also available in three ROX variants – NoROX, LowROX and HighROX.

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Advantages

- Hot-Start technology for better sensitivity
- Accurate genotype calling
- Superior allele clustering
- Compatible with all Real-Time PCR platforms
- Standard and fast cycling conditions

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Applications

- Genotyping single nucleotide polymorphisms (SNPs)
- High throughput genotyping studies

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

Figure

Extracted human genomic DNA was analysed at 2 loci using LNA probes and EliZyme Genotyping MIX. The analysis was performed on a Biorad CFX96 instrument and analysed using the Biorad CFX Manager software version 3.0. Single nucleotide polymorphisms (SNPs) within 2 genes were analysed: IL2RA (A/G) and KYN3 (T/G). Cycling conditions: initial denaturation at 95 °C for 2 minutes, 50 cycles of denaturation at 95 °C for 10 seconds and annealing/extension at 56 °C for 50 seconds. EliZyme Genotyping MIX is fast, accurate and reliable for reproducible allelic discrimination.

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
NEXT GENERATION SEQUENCING

EliZyme Library Quantification Kit

EliZyme Library Quantification Kit for Illumina® contains all the reagents required for successful library quantification. Library quantification is performed by amplification of pre-diluted DNA standards and diluted library samples.

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21

EliZyme Library Quantification Kit

Accurate quantification of library plays a crucial role in next-generation sequencing. When pooling libraries, any deviation of concentration between samples can result in variations in the number of sequenced reads.

Low library concentration may result in low cluster density, whereas high library concentration may result in higher cluster density, which can lead to poor cluster resolution. Both cases result in poor data quality.

EliZyme Library Quantification Kit is based on Real-Time PCR method. EliZyme Library Quantification Kit for Illumina® contains all the reagents required for successful library quantification. Library quantification is performed by amplification of pre-diluted DNA standards and diluted library samples. Primers recognize the adapter sequences used for Illumina® next-generation sequencing. The average Ct value of each DNA standard is plotted against concentration in pM (in log scale) generating a standard curve used for calculating library concentration.

Advantages

- Real-Time PCR for accurate quantification of library
- Accurate and reproducible quantification across a wide range of sample types, concentrations, fragment sizes and GC content
- Non-PCR inhibiting intercalating dye for better signal
- Hot-Start technology

Applications

- Suitable for quantification of NGS libraries for Illumina® platforms containing the p5 and p7 flow cell sequences.
- Compatible with all Real-Time PCR platforms


Figure A
Amplification of EliZyme Library Quantification Kit standards. EliZyme Library Quantification Kit shows equal distribution of amplification curves of the standards compared to the competitor "C" in Fig. A2.


Figure B
Melt curve analysis reveals less primer-dimer formation compared to competitor "C". EliZyme Library Quantification Kit is represented by red curve and competitor "C" by green curve.

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

EliZyme Library Quantification Kit

EZ4815 500 or 1000 rxns (according to reaction vol.)

ACCESSORIES

EliZyme Water for Molecular Biology

EZ9001	2x8 ml
EZ9002	125 ml

EliZyme PCR Water


EZ9101	2x8ml
EZ9102	125 ml

EliZyme Reverse Transcriptase


EZ8005	1x50µl(200U/µl) + 1x200µl buffer	50 rxns
EZ8020	4x50µl (200U/µl) + 4x200µl buffer	200 rxns

EliZyme dNTP MIX

EZ1110	10mM each	1x1 ml
EZ1125	25mM each	1x1 ml

Product portfolio of ELISABETH PHARMACON

Portfolio of EliGene DNA diagnostic kits for Real-Time PCR devices contains products for various applications:

DNA DIAGNOSTICS

- Respiratory Infections
- Sexual Transmitted Infections
- Other Bacterial Infections
- Other Viral Infections

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
- Celiakie
- Spondylitida

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- Colorectum
- GIST
- Lung

EliGene Isolation kits are available in the version for research only or certified for clinical use. In our portfolio there are kits for manual DNA and RNA isolation on column from different materials.

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PORTFOLIO OF KITS

- EliGene MTB Isolation Kit (CE)
- EliGene Urine Isolation Kit (CE)
- EliGene Soil DNA Isolation Kit
- EliGene Blood DNA Isolation Kit
- EliGene Plant DNA Isolation Kit
- EliGene Tissue DNA Isolation Kit
- EliGene Microbial DNA Isolation Kit
- EliGene Viral DNA/RNA Isolation Kit

PORTFOLIO SOUPRAV

- EliGene MTB Isolation Kit (CE)
- EliGene Urine Isolation Kit (CE)
- EliGene Soil DNA Isolation Kit
- EliGene Blood DNA Isolation Kit
- EliGene Plant DNA Isolation Kit
- EliGene Tissue DNA Isolation Kit
- EliGene Microbial DNA Isolation Kit
- EliGene Viral DNA/RNA Isolation Kit

Certificates ELISABETH PHARMACON


ELISABETH PHARMACON has certificated manufacturing and R&D department by ISO 9001:2015 and ISO 13485:2016.

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Elisabeth Pharmacon CZ
Elisabeth Pharmacon UK
Elisabeth Pharmacon SK
members of Elisabeth Pharmacon Group

www.elisabethgroup.com | www.elisabeth.cz | www.eligene.com | www.elizyme.com

Manufacturer:

ELISABETH PHARMACON Ltd.
Rokycanova 4437, 615 00 Brno-Zidenice
phone: +420 542 213 851
e-mail: info@elisabeth.cz
www.elisabeth.cz

Distributor:

EUROPEAN UNION
European Regional Development Fund
Operational Programme Enterprise
and Innovations for Competitiveness